





Perfecto® Mini Random Incidence Performance Data



- E-Mount w/ Fleece Only
- A-Mount w/ 1/2" PET Backing
- A-Mount w/ 2" FBG Backing
- E-Mount w/ 1/2" PET Backing
- A-Mount w/ 1" FBG Backing



| FREQ | E-Mount w/ Fleece Only | E-Mount w/ 1/2" PET Backing | A-Mount w/ 1/2" PET Backing | A-Mount w/ 1" FBG Backing | A-Mount w/ 2" FBG Backing |
|------|---------------------------|--------------------------------|--------------------------------|------------------------------|------------------------------|
| 100 | 0.63 | 0.69 | 0.06 | 0.11 | 0.28 |
| 125 | 0.60 | 0.62 | 0.08 | 0.14 | 0.36 |
| 160 | 0.78 | 0.81 | 0.08 | 0.17 | 0.62 |
| 200 | 0.85 | 0.87 | 0.12 | 0.33 | 0.68 |
| 250 | 0.77 | 0.84 | 0.16 | 0.46 | 0.85 |
| 320 | 0.73 | 0.87 | 0.25 | 0.66 | 1.05 |
| 400 | 0.58 | 0.83 | 0.45 | 0.80 | 1.05 |
| 500 | 0.54 | 0.84 | 0.59 | 1.00 | 1.19 |
| 630 | 0.66 | 0.92 | 0.71 | 1.06 | 1.17 |
| 800 | 0.65 | 0.95 | 0.88 | 1.05 | 1.09 |
| 1000 | 0.66 | 0.97 | 1.01 | 1.03 | 1.08 |
| 1260 | 0.68 | 1.03 | 1.02 | 0.98 | 1.01 |
| 1600 | 0.69 | 1.06 | 0.95 | 0.90 | 0.98 |
| 2000 | 0.69 | 1.02 | 0.80 | 0.84 | 0.91 |
| 2500 | 0.70 | 0.92 | 0.71 | 0.80 | 0.80 |
| 3125 | 0.68 | 0.85 | 0.64 | 0.80 | 0.72 |
| 4000 | 0.64 | 0.79 | 0.60 | 0.78 | 0.64 |
| 5000 | 0.63 | 0.81 | 0.61 | 0.77 | 0.59 |
| NRC | 0.70 | 0.95 | 0.65 | 0.85 | 1.00 |
| SAA | 0.68 | 0.93 | 0.64 | 0.83 | 0.99 |