




Perfecto® Micro MOTIF A Mount Testing Random Incidence Performance Data



— 1" 6LB FBG Backing

— 2" 6LB FBG Backing



| FREQ | 1" 6LB FBG Backing | 2" 6LB FBG Backing |
|------|--------------------|--------------------|
| 100 | 0.17 | 0.41 |
| 125 | 0.29 | 0.65 |
| 160 | 0.34 | 0.91 |
| 200 | 0.51 | 1.14 |
| 250 | 0.70 | 0.90 |
| 320 | 0.87 | 1.13 |
| 400 | 0.98 | 1.06 |
| 500 | 1.09 | 1.05 |
| 630 | 1.11 | 1.00 |
| 800 | 0.93 | 0.91 |
| 1000 | 0.80 | 0.80 |
| 1260 | 0.69 | 0.67 |
| 1600 | 0.58 | 0.55 |
| 2000 | 0.48 | 0.49 |
| 2500 | 0.41 | 0.40 |
| 3125 | 0.32 | 0.34 |
| 4000 | 0.24 | 0.26 |
| 5000 | 0.19 | 0.22 |
| NRC | 0.80 | 0.85 |
| SAA | 0.77 | 0.85 |